

Comparison of in-house PCR and Commercial qPCR for Detecting Gastrointestinal Parasites



Onyimi E.¹, Cunningham L.², Slack I.¹, Jones J¹, Nevin W.D.^{1,3,4}.

1. Clinical Diagnostic Parasitology Lab, Liverpool School of Tropical Medicine, Liverpool UK.
2. Department of Tropical Disease Biology, Liverpool School of Tropical Medicine, Liverpool UK.
3. Academic Department of Military Medicine, Royal Centre of Defence Medicine, Birmingham, UK.
4. Department of Infectious Diseases, Imperial College London, London UK.

Introduction

Gastrointestinal parasite (GIP) infection remains an important global health problem, disproportionally affecting the world's poorest communities. With increased travel and migration there is a drive to improve GIP diagnostics. Microscopic diagnosis requires significant expertise¹, and there is a need for improved data on the validity of molecular diagnostics, particularly for helminths². This study compared faecal microscopy combined with an in-house multiplex parasite qPCR³ to commercially available helminth/protozoa qPCR panels.



Figure 1. Biorad CFX96 PCR Machine

Methods

Concentration microscopy and in-house multiplex qPCR were performed on faecal samples sent to the Liverpool School of Tropical Medicine’s Clinical Diagnostic Parasitology Laboratory (LSTM CDPL) as a reference standard. A selection of faecal samples, both positive and negative for a range of GIPs were then selected to be tested using the commercial Seegene Allplex™ GI-Parasite Assay and/or the commercial Seegene Allplex™ GI-Helminth(I) Assay with the Biorad CFX96 PCR machine. Samples were extracted using the Qiagen DNA extraction kits and stored at -20°C freezers .The in-house multiplex reference standard qPCR panel included *Trichuris*, *Giardia*, *Hookworm spp*, *Strongyloides*, *Schistosoma*, *Ascaris*, *Entamoeba histolytica*, *Entamoeba dispar*, and *Enterobius*.

Results

A total of 125 samples were included. Results from the Seegene™ assays were compared to a reference standard (Microscopy plus in-house qPCR, or microscopy alone in the cases where the GIP target is not included the in-house qPCR panel). For all studied GIP, the Seegene Allplex™ assays showed substantial or near perfect agreement with the reference standard, with Cohen's Kappa ranging from 0.66 for *Ascaris* to 1.00 for *Trichuris*, *Taenia*, *Hymenolepis*, and *Cyclospora*. Sensitivity ranging from 100% for *Trichuris*, *Taenia*, *Hymenolepis*, and *Cyclospora spp* to as low as 50% for *Ascaris*. Specificity was 100% for most of the parasites except *Giardia*, *Strongyloides*, and *Cryptosporidia*. *Hymenolepis*, *Cyclospora* and *Taenia spp* were detected via microscopy alone. The Seegene Allplex™ panels found increased numbers of organisms compared with the reference standard, including four further instances of *Cryptosporidium*, three of *Giardia*, and 24 of *Blastocystis*. There was also 8 instances of microsporidia and nine instances of *Dientamoeba fragilis* not detected by the reference microscopy and in-house qPCR methods. Amongst *Strongyloides* culture positive samples, the seegene panel showed same sensitivity and specificity with in-house qPCR.

Figure 2. Chart showing sensitivity & specificity of each G.I Parasite Sp.




Figure 3. Photograph of the Seegene Allplex™ Result Platform.




Figure 4. Table showing Seegene results against microscopy and in-house PCR reference.

Organism	Ref_pos (n)	Seegene (n)	Sens	Spec	Kappa
Trichuris	6	6	100.00	100.00	1.00
Giardia	11	12	90.91	96.43	0.84
Hookworm	8	7	87.50	100.00	0.93
Strongyloides	15	13	80.00	98.08	0.82
Enterobius	7	6	85.71	100.00	0.92
Ascaris	2	1	50.00	100.00	0.66
Cryptosporidia	8	12	87.50	91.53	0.65
Entamoeba histolytica	11	7	63.64	100.00	0.75
Cyclospora	6	6	100.00	100.00	1.00
Hymenolepis	2	2	100.00	100.00	1.00
Taenia	7	7	100.00	100.00	1.00

Conclusion

The Seegene Allplex™ panels demonstrated good comparability to the reference standard and increased the overall numbers of organisms detected. There was good sensitivity for most helminths (except *Ascaris*, with a small sample size (n=2)). The Seegene Allplex™ overall has a wider panel of organisms and had more identifications compared with the reference in-house methods, including identifying *Dientamoeba fragilis*, *Blastocystis hominis* and microsporidia which are not detected by the in-house qPCR panel. Despite the limitation of a small number of samples for comparison in some cases, the Seegene Allplex qPCR method gave very satisfactory results in this first phase of the comparative analysis. However, a limitation of the Seegene Allplex™ assays is the current absence of *Schistosoma* species, an important global pathogen, from the panels.

References.

1. Buonfrate D, Tamarozzi F, Paradies P, Watts MR, Bradbury RS, Bisoffi Z. The diagnosis of human and companion animal Strongyloides stercoralis infection: challenges and solutions. A scoping review. Advances in parasitology. 2022 Jan 1;118:1-84.
2. Nevin WD, Jones J, Tupper D, Dunbar JA, Wilson D, Ross D, Woolley S, Dodd J, Biswas J, Lamb L, Beeching NJ. Gastrointestinal parasite infections in Nepalese Gurkha recruits arriving in the United Kingdom from 2012–2020. PLOS Neglected Tropical Diseases. 2024 Jan 26;18(1):e0011931.
3. Nevin WD, Cunningham LJ, Mason J, Adams ER, Jones J, Woolley SD, Lamb LE, Beeching NJ, Fletcher TE, O'Shea MK. A comparative study of traditional and molecular diagnostic methods for detection of gastrointestinal parasites in Nepalese migrants to the UK. Journal of Infection. 2024 Dec 1;89(6):106324.